


Cryovial Storage boxes (polycarbonate) 1-2ml Cryotubes



Manufactured from Polycarbonate which is readily able to withstand temperature ranging from 121°C to -196°C makes these storage boxes suitable for various laboratory applications including freezing at ultra-low temperatures.

Features

- Manufactured in Polycarbonate allows ultra-low temperature storage to -196°C
- Designed for storage of Cryovials from 1.0ml to 2.0ml in Liquid Nitrogen or Freezers
- Numeric grid on lid and tray for easy identification of samples
- With vents for storage in liquid nitrogen
- Covers are keyed to the base to prevent misalignment
- 81 well (9 x 9) – **Blue**
- 100 well (10 x 10) – **White**
- Individually shrink-wrapped
- 130mm x 130mm x 53mm (LxWxH)