VACUCARE **Multi-Sample Safety Blood Collection Needles**




VACUCARE **Multi-Sample Safety Blood Collection Needles** are designed to be used with VACUCARE Blood Collection tubes and will fit most other brands.

Needles have triple-bevelled points with bevel up, ensuring smooth venipuncture and venous blood extraction. Hubs and needle guards are colour coded to the needle gauge.

The safety shroud is permanently integrated to the unit. Once deployed the safety shroud is irreversibly locked to and covers the needle, protecting the phlebotomist from needle stick injuries. The locking process can be done single-handedly.

Procedure for locking needle after venipuncture and use



Draw out the needle after venipuncture, perform illustrated actions to lock the needle tube into the safety protective cap

Or by slightly pressing the safety protective cap against the table to lock the needle tube

SPECIFICATION	HUB COLOUR	NEEDLE LENGTH	GAUGE	PACKING
Multi-Sample Safety Blood Collection Needle	Green	38mm or 31mm	21G	48
Multi-Sample Safety Blood Collection Needle	Black	38mm or 31mm	22G	48

5 Portman Square,
London W1H 6NT, UK